



SPECIFICATIONS

Equipment Classification

Class I Equipment, Type B

Room Entry Doorway Requirement

Width: 36 in (91 cm)

Operating Environment (must be STABLE)

Temperature Range: 70°-80°F (21°-27°C)

Relative Humidity: 20%-70% (Non-Condensing)

Altitude: Sea Level to 10,000 ft

Storage Environment

Temperature Range: 40°-100°F (5°-38°C)

BOD POD

Height: 66 in (167 cm)

Width: 33 in (84 cm)

Approx. 45 in (115 cm w/door open)

Depth: 53 in (135 cm)

Approx. 72 in (183 cm w/door open)

Weight: 550 lb (250 kg)

Color: White

Scale

Capacity: 550 lb (250 kg)

Height: 3.5 in (8.9 cm)

Width: 16.0 in (40.6 cm)

Depth: 14.0 in (35.6 cm)

Weight: 23.5 in (10.7 kg)

