

# PRODUCT DATASHEET

## Trolley X1 for COSMED Q-NRG - REF. A-258-900-005

<b>Overall dimensions and weight:</b>	520 mm (W) × 567 mm (D) × 1252 mm (H min) - 12 kg
Work Height (Vesa center)	Fixed (H) 1193 mm
<b>Handle in ABS</b>	Color Blue
<b>Pole with VESA 75-100 attachment</b>	With swivel (-28°/+28°) - 5 Kg max load
<b>Basket</b>	333 mm (int W) × 243 mm (int D) × 101 mm (int H) Transparent PC - 3 Kg max load
<b>Barcode holder</b>	PVC covered steel color White - 1 Kg max load
<b>Aluminium column</b>	
<b>Aluminium base</b>	with Blue bumper
Base footprint	520 mm (W) × 520 mm (D)
4 single wheels with brake	ø 100 mm (4STD)
<b>Optional accessories</b>	
	Easy fix - ref. A-258-500-015
	Single Gas Cylinder Holder - ref. A-258-500-012
<b>Total weight</b>	21 Kg
<b>Packaging dimensions and weight</b>	840 mm (W) × 591 mm (D) × 315 mm (H) - 17 kg
<b>CE Declaration of Conformity</b>	EN 60601-1 Regolamento UE n.2017/745 MDR

